



ADDITIONAL INFORMATION

Weight	1.75 kg
Dimensions	28 × 23 × 14 cm
Carton Qty	8
Pricing Unit	Unit/1
Size	One Size
Min Order Qty	1
ARTG ID	228210
GMDN Code	48047
Medical Device Class	Class 2b
UNSPSC	42172101
Tax Schedule	Non-Taxable
Is Medical Device	Yes
Is Serialized Item	Yes
Is Dangerous Good	Yes

PRODUCT DESCRIPTION

THE HEARTSINE SAMARITAN PAD 360P DEFIBRILLATOR OVERVIEW

- 8-year warranty
- Automatic Shock Delivery
- Easy to Follow Visual and Verbal Guides
- Durable
- Advanced Technology
- Single Battery/Electrode Cartridge

Specially designed for use in [public areas](#), the HeartSine samaritan PAD 360P is a sophisticated fully automatic defibrillator for adult or paediatric use, inside a lightweight and easy-to-operate system. The SAM 360P is a fully automatic defibrillator designed to analyse heart rhythm and automatically deliver an electrical shock (if needed). Its advanced technology is balanced against the demands of real-world use.

How to use the 360P



HEARTSINE SAMARITAN 360P FULLY-AUTOMATIC DEFIBRILLATOR (DG)

PAD-360P

The SAM 360P also uses HeartSine 's practical Pad-Pak™, which incorporates the battery and electrodes in a quick to replace single-use cartridge meaning only one expiration date to monitor and only one item to replace after use. The SAM 360P is suitable for use on adults and children. For children younger than 8 years or weighing less than 25 kg (55 lbs), a Pediatric-Pak™ is available.

- Free replacement Pad-pak after genuine use (conditions apply)
- IP56 Dust and Moisture ingress protection - the highest in the industry
- Metronome for correct speed CPR
- Software updates via USB
- Made in Ireland by the inventors of the first portable **defibrillator**

The HeartSine samaritan PAD 360P is a fully automatic external defibrillator (AED) that analyses the heart rhythm and delivers an electrical shock to victims of Sudden Cardiac Arrest in order to restore the heart to normal rhythm. This user-friendly AED provides easy to follow visual and audio prompts to guide the rescuer through pad application. Shock delivery, if required, is fully automatic which means there is no shock button to press.

[Register your Product Warranty Here.](#)

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Defibrillator Samaritan PAD 360P